


CATGUT

CHROM



Type	Twisted multifilament with a monofilament appearance
Composition	Strands of purified collagen taken from the serosal layer of selected bovines, tanned with chrome salts. Origin from BSE-free classified countries.
Coating	None
Colour	Dark brown
Tissue reaction	Moderate
Absorption	By phagocytosis. The rate of absorption is slowed down by chromizing. Complete mass absorption in approximately 90 days. Approximately 50% of tensile strength remain after 14 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

CATGUT

CHROM

CHARACTERISTICS

- ✓ High tensile strength
- ✓ Absorption by phagocytosis
- ✓ Twisted structure
- ✓ Uniform diameter by centerless polishing (monofilament appearance)
- ✓ Hermetically sealed packing

BENEFITS

Excellent tissue support

Uniform absorption

Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

Soft passage through the tissues without sawing, tissue drag and trauma
Excellent knot tie down

Guaranteed seal and product sterility

INDICATIONS

General soft tissue closing and/or ligation

especially in

- General surgery
- Gastrointestinal surgery
- Gynaecology
- Obstetrics
- Urology
- Ophthalmic surgery

CONTRAINDICATIONS

Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).


In the region of infected wounds

Known sensitivities or allergies to collagen or chromium

TENSILE STRENGTH



ABSORPTION





Without needle

Metric USP Length Order number

2	4/0	150 cm	<i>brown</i>	2200015
3	3/0	150 cm	<i>brown</i>	2300015
3.5	2/0	150 cm	<i>brown</i>	2350015
4	0	150 cm	<i>brown</i>	2400015
5	1	150 cm	<i>brown</i>	2500015
6	2	150 cm	<i>brown</i>	2600015
7	3	150 cm	<i>brown</i>	2700015
8	4	150 cm	<i>brown</i>	2800015
10	6	150 cm	<i>brown</i>	21000015



Round bodied taper point

Metric USP Length Order number

1/2 circle



HR 15 mm



1	6/0	75 cm	<i>brown</i>	2100115
1.5	5/0	75 cm	<i>brown</i>	2150115
2	4/0	75 cm	<i>brown</i>	2200115

HR 20 mm



1.5	5/0	75 cm	<i>brown</i>	2150120
2	4/0	75 cm	<i>brown</i>	2200120
3	3/0	75 cm	<i>brown</i>	2300120
3.5	2/0	75 cm	<i>brown</i>	2350120
4	0	75 cm	<i>brown</i>	2400120

HR 26 mm



1.5	5/0	75 cm	<i>brown</i>	2150126
2	4/0	75 cm	<i>brown</i>	2200126
3	3/0	75 cm	<i>brown</i>	2300126
3.5	2/0	75 cm	<i>brown</i>	2350126
4	0	75 cm	<i>brown</i>	2400126
5	1	75 cm	<i>brown</i>	2500126

HR 30 mm







2	4/0	75 cm	<i>brown</i>	2200130
3	3/0	75 cm	<i>brown</i>	2300130
3.5	2/0	75 cm	<i>brown</i>	2350130
4	0	75 cm	<i>brown</i>	2400130
5	1	75 cm	<i>brown</i>	2500130



HR 37 mm















2	4/0	75 cm	<i>brown</i>	2200137
3	3/0	75 cm	<i>brown</i>	2300137
3.5	2/0	75 cm	<i>brown</i>	2350137
3.5	2/0	90 cm	<i>brown</i>	2350137-90
4	0	75 cm	<i>brown</i>	2400137
5	1	75 cm	<i>brown</i>	2500137

 Round bodied taper point		Metric	USP	Length		Order number
 1/2 circle	HR 40 mm	3.5	2/0	75 cm	<i>brown</i>	2350140
		3.5	2/0	90 cm	<i>brown</i>	2350140-90
		4	0	75 cm	<i>brown</i>	2400140
		5	1	75 cm	<i>brown</i>	2500140
		5	1	90 cm	<i>brown</i>	2500140-90
		6	2	75 cm	<i>brown</i>	2600140
	HR 45 mm	4	0	75 cm	<i>brown</i>	2400145
		6	2	75 cm	<i>brown</i>	2600145
	HR 48 mm	5	1	75 cm	<i>brown</i>	2500148
	HR 50 mm	6	2	75 cm	<i>brown</i>	2600150
		6	2	90 cm	<i>brown</i>	2600150-90

 Round bodied taper point reverse cutting double armed		Metric	USP	Length		Order number
 1/2 circle	HR 37 mm HS 37 mm	3.5	2/0	150 cm	<i>brown</i>	2350137-1437

 Blunt		Metric	USP	Length		Order number
 3/8 circle	DRB 63 mm	5	1	75 cm	<i>brown</i>	2501963

 Reverse cutting		Metric	USP	Length		Order number
3/8 circle 	DS 12 mm		1	6/0	45 cm	<i>brown</i> 2101512
			1.5	5/0	45 cm	<i>brown</i> 2151512
	DS 16 mm		1.5	5/0	75 cm	<i>brown</i> 2151516
			2	4/0	75 cm	<i>brown</i> 2201516
	DS 19 mm		1.5	5/0	75 cm	<i>brown</i> 2151519
			2	4/0	75 cm	<i>brown</i> 2201519
			3	3/0	75 cm	<i>brown</i> 2301519
	DS 24 mm		3.5	2/0	75 cm	<i>brown</i> 2351519
			2	4/0	75 cm	<i>brown</i> 2201524
			3	3/0	75 cm	<i>brown</i> 2301524
			3.5	2/0	75 cm	<i>brown</i> 2351524
			4	0	75 cm	<i>brown</i> 2401524
	DS 30 mm		5	1	75 cm	<i>brown</i> 2501524
			3	3/0	75 cm	<i>brown</i> 2301530
			3.5	2/0	75 cm	<i>brown</i> 2351530
			4	0	75 cm	<i>brown</i> 2401530
	DS 37 mm		5	1	75 cm	<i>brown</i> 2501530
			3	3/0	75 cm	<i>brown</i> 2301537
			3.5	2/0	75 cm	<i>brown</i> 2351537
			4	0	75 cm	<i>brown</i> 2401537
	DS 40 mm		5	1	75 cm	<i>brown</i> 2501537
			3.5	2/0	75 cm	<i>brown</i> 2351540
			6	2	75 cm	<i>brown</i> 2601540
			7	3	75 cm	<i>brown</i> 2701540

 Reverse cutting		Metric	USP	Length		Order number
1/2 circle 	HS 37 mm		4	0	75 cm	<i>brown</i> 2401437
			5	1	75 cm	<i>brown</i> 2501437