

CATGUT

PLAIN



Type	Twisted multifilament with a monofilament appearance
Composition	Strands of purified collagen taken from the serosal layer of selected bovines. Origin from BSE-free classified countries.
Coating	None
Colour	Ivory
Tissue reaction	Moderate
Absorption	By phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

CATGUT

PLAIN

CHARACTERISTICS

- ✓ High tensile strength
- ✓ Absorption by phagocytosis
- ✓ Twisted structure
- ✓ Uniform diameter by centerless polishing (monofilament appearance)
- ✓ Hermetically sealed packing

BENEFITS

Excellent tissue support

Uniform absorption

Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

Soft passage through the tissues without sawing, tissue drag and trauma
Excellent knot tie down

Guaranteed seal and product sterility

INDICATIONS

General soft tissue closing and/or ligation

especially in

- General surgery
- Gastrointestinal surgery
- Gynaecology
- Obstetrics
- Urology
- Ophthalmic surgery

CONTRAINDICATIONS

Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).


In the region of infected wounds

Known sensitivities or allergies to collagen

TENSILE STRENGTH



ABSORPTION





Without needle

Metric USP Length Order number

2	4/0	150 cm	ivory	1200015
3	3/0	150 cm	ivory	1300015
3.5	2/0	150 cm	ivory	1350015
4	0	150 cm	ivory	1400015
5	1	150 cm	ivory	1500015
6	2	150 cm	ivory	1600015
7	3	150 cm	ivory	1700015











Round bodied taper point

Metric USP Length Order number

1/2 circle



HR 16 mm		1.5	5/0	75 cm	ivory	1150116
		2	4/0	75 cm	ivory	1200116
HR 20 mm		2	4/0	75 cm	ivory	1200120
		3	3/0	75 cm	ivory	1300120
		3.5	2/0	75 cm	ivory	1350120
HR 26 mm		2	4/0	75 cm	ivory	1200126
		3	3/0	75 cm	ivory	1300126
		3.5	2/0	75 cm	ivory	1350126
		4	0	75 cm	ivory	1400126
HR 30 mm		3	3/0	75 cm	ivory	1300130
HR 37 mm		3	3/0	75 cm	ivory	1300137
		3.5	2/0	75 cm	ivory	1350137
		4	0	75 cm	ivory	1400137
		4	0	90 cm	ivory	1400137-90
		5	1	75 cm	ivory	1500137
HR 40 mm		5	1	90 cm	ivory	1500140-90
HR 45 mm		4	0	75 cm	ivory	1400145
		5	1	75 cm	ivory	1500145
		6	2	75 cm	ivory	1600145
HR 50 mm		3.5	2/0	75 cm	ivory	1350150



Reverse cutting

Metric






USP

Length

Order number

3/8 circle



DS 12 mm		1.5	5/0	45 cm	<i>ivory</i>	1151512
DS 16 mm		1	6/0	45 cm	<i>ivory</i>	1101516
		1.5	5/0	75 cm	<i>ivory</i>	1151516
		2	4/0	75 cm	<i>ivory</i>	1201516
		3	3/0	75 cm	<i>ivory</i>	1301516
DS 19 mm		2	4/0	75 cm	<i>ivory</i>	1201519
		3	3/0	75 cm	<i>ivory</i>	1301519
		3.5	2/0	75 cm	<i>ivory</i>	1351519
DS 24 mm		2	4/0	75 cm	<i>ivory</i>	1201524
		3	3/0	75 cm	<i>ivory</i>	1301524
		3.5	2/0	75 cm	<i>ivory</i>	1351524
		4	0	75 cm	<i>ivory</i>	1401524
DS 37 mm		3.5	2/0	75 cm	<i>ivory</i>	1351537